TACTICAL FIELD TRAUMA KIT W/QUIKCLOT - TACTICAL FIELD TRAUMA KIT WITH QUIKCLOT

Nobody thinks that they'll be the ones to suffer life-threatening injuries when they go into the field, so oftentimes, many hikers, hunters and outdoorsmen leave essential medical supplies at their base camp or in their vehicle. The Adventure Medical Tactical Field Trauma Kit with QuikClot ensures that these adventurers always have essential medical equipment with them when disaster strikes. The Adventure Medical Tactical Field Trauma Kit with QuikClot is sized and weighted right to be easily added to any hiking pack while containing all of the gear necessary to handle everything from basic cuts all the way to enhanced wound care. Included QuikClot Sport dressing stops heavy bleeding fast Designed for use by up to 4 people Provides supplies for 1- to 2-day treks The Adventure Medical Tactical Field Trauma Kit with QuikClot also includes all of the materials necessary to stabilize hurt individuals and keep them out of pain until help arrives.

Attributes

- Name: TACTICAL FIELD TRAUMA KIT WITH QUIKCLOT
- Manufacturer: ADVENTURE MEDICAL KITS
- Product no.: 100025344
- Mfr. No.: 2064-0291
- Size: 4 Person
- Delivery weight: 0.363kg
- UPC: 707708102912

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [English: TACTICAL FIELD TRAUMA KIT W/QUIKCLOT Safety Instruction Guide](#)
- [Svenska: Säkerhetsinstruktioner för TACTICAL FIELD TRAUMA KIT med QuikClot](#)
- [Český: Taktická traumatická souprava Adventure Medical s QuikClot Bezpečnostní pokyny pro uživatele](#)

TACTICAL FIELD TRAUMA KIT W/QUIKCLOT Safety Instruction Guide

Introduction

Thank you for choosing the Adventure Medical Tactical Field Trauma Kit with QuikClot. This kit is designed to provide essential medical supplies for outdoor adventurers, ensuring safety during unexpected emergencies. This guide outlines important safety instructions and guidelines to help you use the kit effectively and safely.

General Safety Guidelines

- Ensure that all users of the kit are familiar with its contents and how to use each item.
- Store the kit in a cool, dry place, away from direct sunlight and moisture.
- Regularly check the kit for expired items and replace them as necessary.
- Keep the kit out of reach of children when not in use.
- In case of an emergency, always call for professional medical assistance as soon as possible.

Specific Safety Precautions for Use

- Use the QuikClot dressing only as directed to control bleeding. Improper use may lead to further injury.
- Avoid using the kit on individuals with known allergies to any of the materials included.
- Do not use the kit as a substitute for professional medical care.
- If you are unsure about how to use any item in the kit, seek guidance from a trained professional.
- Always wash your hands before and after treating any injuries to prevent infection.

Instructions for Installation and Usage

Storage

- Keep the kit in an accessible location in your hiking pack or vehicle.
- Ensure that the kit is securely closed to prevent loss of contents.

Using the Kit

1. **Assess the Situation:** Check the injured person's condition and call for emergency help if necessary.
2. **Control Bleeding:**
 - If there is heavy bleeding, apply the QuikClot dressing directly to the wound.
 - Apply firm pressure for a minimum of three minutes.
3. **Stabilize Injuries:**
 - Use splints or bandages included in the kit to stabilize any fractures or sprains.
 - Ensure that the injured area is immobilized until professional help arrives.
4. **Monitor the Injured Person:**
 - Keep the individual calm and comfortable.
 - Monitor their vital signs until help arrives.

PostUse Care

- After use, dispose of any singleuse items according to local regulations.
- Replenish any used items in the kit to ensure it is ready for future use.

Disposal Instructions

- Dispose of used items safely and in accordance with local waste management regulations.
- Do not throw away items in regular trash if they are biohazardous or contain bodily fluids.

Contact Information for Further Support

For any questions or concerns regarding the Adventure Medical Tactical Field Trauma Kit with QuikClot, please contact the EUs-based support team. Always refer to the product packaging or official website for the most current information.

Conclusion

The Adventure Medical Tactical Field Trauma Kit with QuikClot is an essential tool for anyone venturing into the outdoors. By following these safety instructions, you can ensure that you are prepared for emergencies while minimizing risks. Always prioritize safety and seek professional medical assistance when needed. Thank you for your commitment to safety during your outdoor adventures.

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