HEMOSTATIC GRANULES, 15G PACK

Celox™ Granules can stop severe bleeding fast. The Celox™ Granules mold to the shape of the wound to get firm pressure on to the source of the bleeding even in complex shaped injuries. Safe, easy to use, cost effective and capable of applying to a variety of different wounds No Exothermic reaction – Does not generate heat Works in hypothermic conditions and clots heparinized blood Simple operation – Tear open, pour, press, pack Celox™ has been proven to work on moderate to severe bleeding and has been extensively tested. Celox granules are available in a 15g sachet for situations where every gram of weight saved is vital, still with enough material to treat a severe injury.

Attributes

- Name: HEMOSTATIC GRANULES, 15G PACK
- Manufacturer: CELOX MEDICAL
- Product no.: 100025879
- Mfr. No.: V12090
- Size: 1 Person
- Delivery weight: 0.005kg
- UPC: 755702981620

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Safety Instruction Guide for HEMOSTATIC GRANULES

Introduction

HEMOSTATIC GRANULES, specifically Celox™ Granules, are designed to stop severe bleeding quickly and effectively. This guide provides essential safety instructions to ensure safe use, compliance with EU General Product Safety Regulation (GPSR), and optimal effectiveness of the product.

General Safety Guidelines

- Always read and understand the instructions before use.
- Keep out of reach of children to prevent accidental ingestion.
- Use only as directed. This product is intended for external use only.
- Ensure the product is stored in a cool, dry place away from direct sunlight.
- Check the product for any damage or tampering before use.
- If you experience any adverse reactions, discontinue use and seek medical attention.

Specific Safety Precautions for Use

- Do not use Celox™ Granules on deep puncture wounds or animal bites.
- Avoid contact with eyes and mucous membranes. If contact occurs, rinse thoroughly with water.
- Do not use on patients with known allergies to chitosan or other components of the granules.
- In case of severe bleeding, if the bleeding does not stop after applying the granules, seek medical assistance immediately.
- Do not use in combination with other hemostatic agents unless directed by a healthcare professional.

Instructions for Installation and Usage

1. Preparation:

- Ensure the wound is clean. Remove any debris or foreign objects if possible.
- Wear gloves if available to maintain hygiene.

2. Application:

- Tear open the sachet of Celox™ Granules.
- Pour the granules directly onto the wound, covering the entire area.
- Press firmly on the granules to ensure they mold to the shape of the wound.
- Maintain pressure on the wound for at least 23 minutes to facilitate clotting.

3. PostApplication:

- Once bleeding has stopped, cover the wound with a sterile dressing or bandage.
- Dispose of any used materials properly.

Disposal Instructions

- Dispose of any unused or expired Celox™ Granules in accordance with local regulations.
- Do not dispose of the product in waterways or general waste unless specified by local waste management guidelines.

Contact Information for Further Support

For any safety inquiries or further information regarding the use of HEMOSTATIC GRANULES, please refer to the product packaging or manufacturer's website for contact details.

By following these guidelines, you can ensure the safe and effective use of HEMOSTATIC GRANULES. Always prioritize safety and consult a healthcare professional if you have any concerns regarding the use of this product.

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