RAPID HEMOSTATIC Z-FOLD GAUZE - 3"X5' CELOX RAPID

HEMOSTATIC Z-FOLD GAUZE

Time is critical when treating severe injuries in the field and saving time increases the chance of survival for the casualty in an emergency or hostile situation. To be effective in real use, hemostatic dressings need to work fast. The low profile pack means Celox Rapid Ribbon is easily concealed and the reduced size makes it ideal for treating narrow entry wounds such as bullet, blast fragment and stab injuries. Works with just 60 seconds compression. Stops severe bleeding – fast. Rapid action reduces blood loss. Stops hypothermic bleeding as found in severe traumatic injuries. The current generation of hemostatic agents need at least 3 minutes of compression to work. Celox Rapid Ribbon Gauze removes the delay by speeding up packing time and reducing compression time to just 60 seconds while using the proven hemostatic technology of Celox.

Attributes

- Name: 3"X5' CELOX RAPID HEMOSTATIC Z-FOLD GAUZE
- Manufacturer: CELOX MEDICAL
- Product no.: 100025884
- Mfr. No.: TS-ZR-035
- Size: 1 Person
- Delivery weight: 0.005kg
- UPC: 755702981835

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Safety Instruction Guide for RAPID HEMOSTATIC ZFOLD GAUZE 3"X5' CELOX

Introduction

The RAPID HEMOSTATIC ZFOLD GAUZE 3"X5' CELOX is designed for immediate use in emergency situations to control severe bleeding. This guide provides essential safety instructions to ensure effective and safe use of the product in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Always use the hemostatic gauze as directed to minimize risks.
- Ensure that the product is stored in a cool, dry place away from direct sunlight.
- Regularly check the packaging for any signs of damage or expiration.
- Keep out of reach of children and ensure that it is used by trained personnel only.
- In case of any adverse reactions or complications, seek medical attention immediately.

Specific Safety Precautions for Use

- Use the gauze only for its intended purpose: to control severe bleeding.
- Do not use if the packaging is damaged or if the product has expired.
- Ensure that the area of application is clean to avoid infection.
- Apply direct pressure to the wound after placing the gauze to ensure effective hemostasis.
- If bleeding does not stop within 60 seconds of application, seek immediate medical help.
- Be aware of potential allergic reactions; check for any known allergies to the components of the gauze.

Instructions for Installation and Usage

1. Preparation:

- Wash your hands thoroughly before handling the gauze.
- Ensure that the area around the wound is as clean as possible.

2. Application:

- Remove the gauze from its packaging carefully.
- Place the gauze directly on the wound.
- Apply firm pressure for 60 seconds to allow the hemostatic agents to work effectively.
- If necessary, secure the gauze with additional bandaging to maintain pressure.

3. PostApplication:

- Monitor the wound for any signs of continued bleeding.
- If bleeding persists after 60 seconds, do not remove the gauze; seek medical assistance immediately.
- Dispose of used gauze in accordance with local regulations regarding medical waste.

Disposal Instructions

- Dispose of any used or expired gauze in accordance with local medical waste disposal regulations.
- Do not reuse the gauze once it has been used for bleeding control.

Contact Information for Further Support

For any safety inquiries or to report unsafe products, please refer to the appropriate EU contact points for consumer safety.

This safety instruction guide is designed to provide clear and concise information for the effective use of RAPID HEMOSTATIC ZFOLD GAUZE 3"X5' CELOX. Always prioritize safety and seek professional assistance when needed.

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