# FOXSEAL CHEST SEAL OCCLUSIVE DRESSING 2/PACK

Foxseal Chest Seal is a pioneering occlusive dressing for the treatment of open chest wounds as a tool in prevention of pneumothorax. Compact, fits IFAK without folding, to avoid creases. Strong adhesion in extreme conditions. Two per pack for ballistic wounds. Simple no-valve design. Vent with large, easy to operate tab. Other dressing options are often folded when carried in kits, leading to a permanent crease which can prevent a true seal. Foxseal is a compact size to fit neatly into pouches, to avoid creasing. Two seals per pack allow entry and exit wounds to be treated. Patient monitoring is vital in these injury types and the design is a simple valve-free option, for simple use. The chest can be vented with the easily-gripped tab and the adhesive re-seals.

## Attributes

- Name: FOXSEAL CHEST SEAL OCCLUSIVE DRESSING 2/PACK
- Manufacturer: CELOX MEDICAL
- Product no.: 100025890
- Mfr. No.: TS-FSOD1
- Size: 2 Person
- Delivery weight: 0.005kg
- UPC: 755702981873

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# FOXSEAL CHEST SEAL OCCLUSIVE DRESSING

## SAFETY INSTRUCTION GUIDE

### Introduction

Thank you for choosing the FOXSEAL Chest Seal Occlusive Dressing. This product is designed for the effective treatment of open chest wounds and to help prevent pneumothorax. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

### General Safety Guidelines

- Ensure the product is used only for its intended purpose.
- Keep out of reach of children to prevent accidental use.
- Check the packaging for any signs of damage or tampering before use.
- Follow all instructions carefully to avoid potential hazards.
- Report any adverse reactions or incidents to the appropriate authorities.

### Specific Safety Precautions for Use

- Verify that the dressing is sealed and undamaged before use.
- Do not use the dressing if it has passed its expiration date.
- Avoid contact with open wounds if you have not washed your hands or are not wearing gloves.
- Monitor the patient closely for signs of respiratory distress or other complications.
- Use the dressing only in accordance with established first aid procedures.

### Instructions for Installation and Usage

#### 1. Preparation:

- Ensure the area is clean and dry.
- Wash hands thoroughly or wear disposable gloves.

#### 2. Application:

- Remove the FOXSEAL Chest Seal from its packaging.
- Peel off the backing to expose the adhesive side.
- Place the dressing over the wound, ensuring a tight seal around the edges.
- Ensure that the vent tab is positioned upward and is easily accessible.

#### 3. Ventilation:

- If the patient shows signs of respiratory distress, use the vent tab to release trapped air.
- Gently pull the tab to allow air to escape while maintaining the seal around the wound.

#### 4. Monitoring:

- Continuously monitor the patient's condition for any changes.
- If the dressing becomes loose or the patient experiences worsening symptoms, replace the dressing as necessary.

#### 5. Removal:

- Carefully peel the dressing off the skin, avoiding any further injury.
- Dispose of the used dressing according to local regulations.

### Disposal Instructions

- Dispose of the FOXSEAL Chest Seal in accordance with local waste management regulations.
- Do not reuse the dressing.
- Ensure that used dressings are placed in a sealed bag to prevent contamination.

## **Contact Information for Further Support**

For safety inquiries or further assistance, please refer to the manufacturer's contact information provided on the packaging or product insert.

By following these guidelines, you can ensure the safe and effective use of the FOXSEAL Chest Seal Occlusive Dressing. Thank you for prioritizing safety and health in your first aid practices.

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