# RIGID TQ CASE FOR C-A-T WITH TEK-LOK - RIGID CAT GEN 7 CASE W/TEK-LOK, MULTICAM

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.

## Attributes

- Name: RIGID CAT GEN 7 CASE W/TEK-LOK, MULTICAM
- Manufacturer: ELEVEN 10
- Product no.: 100043715
- Mfr. No.: E103021BMTC
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615672028

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# Rigid TQ Case for CAT with TekLok Safety Instructions

## Introduction

Thank you for purchasing the Eleven 10® Rigid TQ Case® for the CAT® Gen 7. This hardsided tourniquet case is designed for quick access to your tourniquet in emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

## General Safety Guidelines

The following guidelines are designed to ensure your safety and the proper use of the Rigid TQ Case:

- Always inspect the case and tourniquet before use to ensure they are in good condition.
- Keep the case clean and free from debris to ensure proper function.
- Store the case in a location that is easily accessible in an emergency.
- Ensure that the case is securely mounted to your duty belt, vest, or holster to prevent accidental detachment.
- Follow all local regulations and guidelines regarding the use of tourniquets and first aid supplies.

## Specific Safety Precautions for Use

When using the Rigid TQ Case, please observe the following precautions:

- Ensure that the tourniquet is properly seated in the case and secured using friction.
- Avoid using the case if it has visible damage or wear that may compromise its integrity.
- Do not attempt to modify or alter the case in any way.
- Use the case only for its intended purpose, which is to hold and deploy a CAT® Tourniquet.
- Keep the case out of reach of children and vulnerable individuals to prevent misuse.

## Instructions for Installation and Usage

Follow these steps to install and use the Rigid TQ Case effectively:

### 1. Mounting the Case:

- Choose a suitable mounting location (duty belt, MOLLE vest, or Safariland® drop leg holster).
- If mounting on a duty belt, ensure the BladeTech TEKLOK Belt Attachment is securely fitted to belts up to 2.25 inches wide.
- For MOLLE/PALS compatible vests, attach the case using compatible mounting solutions.
- Ensure the case is mounted at a height that allows for easy access.

### 2. Inserting the Tourniquet:

- Open the case and insert the CAT® Tourniquet, ensuring it is seated properly.
- Close the case securely, ensuring the tourniquet is held in place by friction.

### 3. Deploying the Tourniquet:

- In an emergency, locate the exposed portion of the tourniquet.
- Pull swiftly to release the tourniquet from the case.
- Follow the tourniquet's instructions for application.

### 4. Regular Maintenance:

- Periodically check the case and tourniquet for any signs of wear or damage.
- Clean the case with a damp cloth and mild detergent as needed.

## Disposal Instructions

When the Rigid TQ Case or the tourniquet reaches the end of its life cycle or is no longer usable:

- Dispose of the case and tourniquet in accordance with local regulations for hazardous waste.
- Do not throw the product in regular trash if it contains materials that may be considered hazardous.

## **Contact Information for Further Support**

For any inquiries regarding safety, usage, or product issues, please contact the manufacturer or your local distributor. Ensure you have the product information available for reference.

Please remember to report any unsafe products or accidents to the relevant authorities and check for recall updates on the EU's Safety Gate platform. Your safety is our priority.

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