

RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, BLACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, BLACK](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043721
- Mfr. No.: E10-3011B-BLK
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668236

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [RIGID TQ Case Safety Instruction Guide](#)
- [About Us](#)

RIGID TQ Case Safety Instruction Guide

Introduction

Thank you for choosing the RIGID TQ Case for SOF®TT/SOF®TTW by Eleven 10®. This guide provides important safety instructions and usage information to ensure the safe and effective use of your tourniquet case. Please read this document carefully before using the product.

General Safety Guidelines

- Ensure the product is used only for its intended purpose, which is to securely hold and provide rapid access to SOF®TT or SOF®TTWide tourniquets.
- Regularly inspect the RIGID TQ Case for any signs of damage or wear. If any damage is found, discontinue use immediately.
- Keep the product out of reach of children to prevent accidental use or injury.
- Be aware of your surroundings when deploying the tourniquet to avoid injury to yourself or others.
- Follow all local laws and regulations regarding the carrying and use of tourniquets.

Specific Safety Precautions for Use

- Always ensure the tourniquet is properly secured in the case before use.
- Use the case only with compatible tourniquets (SOF®TT or SOF®TTWide).
- Do not attempt to modify or alter the RIGID TQ Case in any way.
- Avoid exposing the case to extreme temperatures or harsh chemicals, as this may compromise its integrity.
- When deploying the tourniquet, use onehanded gross motor skills as designed to ensure quick access, especially in emergencies.

Instructions for Installation and Usage

Installation

- The RIGID TQ Case can be mounted in various positions. Choose the mounting option that best suits your needs:
 - Police Duty Belt (using BladeTech TEKLOK, fits up to a 2.25" belt)
 - MOLLE/PALS Attachment (using Tactical Tailor Malice Clip)
 - Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately)
- Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained. Instructions are included with the product.

Usage

- To deploy the tourniquet, grab the exposed portion firmly and pull swiftly to release it from the case.
- Practice deploying the tourniquet in a safe environment to familiarize yourself with the process.
- Ensure that the tourniquet is positioned correctly on the limb requiring treatment as per the manufacturer's instructions.

Disposal Instructions

- Dispose of the RIGID TQ Case in accordance with local regulations for plastic products.
- If the case is damaged beyond repair, ensure it is disposed of in a way that prevents any potential misuse.

Contact Information for Further Support

For any questions or concerns regarding the RIGID TQ Case, please refer to the product packaging or visit the manufacturer's website for additional support and resources.

This safety instruction guide is designed to help you use the RIGID TQ Case safely and effectively. Always prioritize safety and consult professionals if you are unsure about any aspect of the product's use. Thank you for your attention, and stay safe!

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk