

# RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, ORANGE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



## Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, ORANGE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043728
- Mfr. No.: E10-3011M-ORG
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668397

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Rigid TQ Case Safety Instructions](#)
- [About Us](#)

# Rigid TQ Case Safety Instructions

## Introduction

Thank you for choosing the Rigid TQ Case for the SOF®TT/SOF®TTW. This safety instruction guide is designed to help you understand the proper use, installation, and disposal of your tourniquet case. Please read these instructions carefully to ensure safe and effective use of the product.

## General Safety Guidelines

To maximize your safety and the effectiveness of the Rigid TQ Case, please adhere to the following guidelines:

- Ensure that the case is used only for its intended purpose: to hold SOF®TT or SOF®TTW tourniquets.
- Regularly inspect the case for any signs of damage or wear. If damaged, discontinue use immediately.
- Keep the case out of reach of children and vulnerable individuals to prevent accidental use.
- Follow all local regulations regarding the use and carrying of medical devices.
- In case of an emergency, always prioritize calling for professional medical help.

## Specific Safety Precautions for Use

When using the Rigid TQ Case, please observe the following safety precautions:

- Ensure that the tourniquet is securely held in place within the case before use.
- Practice deploying the tourniquet from the case in a safe environment to familiarize yourself with the process.
- Avoid using excessive force when pulling the tourniquet from the case to prevent damage.
- If the case becomes wet, dry it thoroughly before storing or using it again.
- Do not modify or alter the case in any way, as this may compromise its safety and effectiveness.

## Instructions for Installation and Usage

Follow these steps for proper installation and usage of the Rigid TQ Case:

### 1. Mounting the Case:

- Choose a suitable mounting location, such as a duty belt or MOLLE/PALS compatible vest.
- For police duty belt installation, ensure the BladeTech TEKLOK Belt Attachment fits securely onto a belt up to 2.25 inches wide.
- For MOLLE/PALS installation, attach the Tactical Tailor Malice Clip to the designated area on your vest.

### 2. Placing the Tourniquet:

- Insert the SOF®TT or SOF®TTW tourniquet into the case, ensuring it is seated properly.
- Make sure the exposed portion of the tourniquet is easily accessible for rapid deployment.

### 3. Using the Tourniquet:

- In an emergency, quickly grab the exposed portion of the tourniquet.
- With a swift pull, release the tourniquet from the case.
- Follow the manufacturer's instructions for applying the tourniquet effectively.

## **Disposal Instructions**

When the Rigid TQ Case reaches the end of its life cycle or if it is damaged beyond repair, please dispose of it responsibly:

- Check local regulations for proper disposal methods of plastic products.
- If the case is recyclable, separate it from other waste and place it in the appropriate recycling bin.
- If not recyclable, dispose of the case in accordance with local waste management guidelines.

## **Contact Information for Further Support**

For any safety inquiries or further assistance regarding the Rigid TQ Case, please refer to the contact information provided with your product packaging or visit the manufacturer's website.

Thank you for prioritizing safety and compliance. Your attention to these guidelines will help ensure effective use of the Rigid TQ Case in critical situations.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)