MICRO TRAUMA KIT NOW! ADVANCED SUPPLIES

The Blue Force Gear Micro Trauma Kit Advanced refill is not just another small tactical first aid kit. Each component was specifically chosen for its combat proven effectiveness combined with never before seen packaging techniques to make them as small as possible. These supplies are everything necessary to address the “Big Three” when it comes to preventable deaths due to traumatic injury; hemorrhage, occluded airway, and tension pneumothorax. The Micro Trauma Kit Advanced supplies were designed from the ground up around a specific need for an ultra-concealable trauma kit. The Kit contains the following: 1x CLEER Medical Trauma Bandage (4" x 24" Flat Packed compression dressing) 1x Celox Rapid (3" x 24" Z-Fold) 2x FOXSeal Occlusive Dressing 1x North American Rescue Tan Bear Claw Gloves 4x CLEER Medical Easy Tape 1x ARS for Needle Decompression (14 gauge x 3.25") 1x Rusch Nasopharyngeal Airway (28FR with lubricating jelly) ***Brands and/or packaging of supplies may vary from items in photos***

Attributes

- Name: MICRO TRAUMA KIT NOW! ADVANCED SUPPLIES
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100044829
- Mfr. No.: P-TKR-03
- Delivery weight: 0.2kg
- UPC: 812114025745

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Safety Instruction Guide for MICRO TRAUMA KIT NOW! ADVANCED SUPPLIES

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! ADVANCED SUPPLIES. This kit is designed to provide essential medical supplies for addressing traumatic injuries effectively. It is crucial to understand how to use this kit safely and effectively. This guide outlines important safety instructions and precautions to ensure your safety and the safety of others.

General Safety Guidelines

- Always read and understand the instructions before using the kit.
- Keep the kit out of reach of children and vulnerable individuals.
- Inspect the kit regularly to ensure all components are present and in good condition.
- Store the kit in a cool, dry place, away from direct sunlight.
- Do not use expired items; check expiration dates regularly.
- Report any incidents or unsafe products to the appropriate authorities.

Specific Safety Precautions for Use

- **Hemorrhage Control:** Use the CLEER Medical Trauma Bandage and Celox Rapid as directed to control bleeding. Ensure you apply adequate pressure and follow the instructions for use.
- **Airway Management:** Use the Rusch Nasopharyngeal Airway only if you are trained in airway management. Improper use can cause injury.
- **Tension Pneumothorax:** The ARS for Needle Decompression should only be used by trained personnel. Misuse can lead to serious complications.
- **Personal Protective Equipment:** Always wear gloves (included in the kit) when providing care to prevent infection.
- **Occlusive Dressing:** Use FOXSeal Occlusive Dressing to seal open wounds and prevent air entry. Ensure it is applied correctly to avoid further complications.

Instructions for Installation and Usage

1. Preparation:

- Ensure you have a clear and safe area to provide care.
- Wash your hands thoroughly before handling any supplies.
- Gather necessary items from the kit based on the injury type.

2. Using the Kit:

- **For Bleeding:**
 - Open the CLEER Medical Trauma Bandage and apply it directly to the wound.
 - Apply pressure and secure the bandage in place.
 - If bleeding does not stop, consider applying additional dressings and seek medical help.
- **For Airway Issues:**
 - If trained, gently insert the Rusch Nasopharyngeal Airway into the nasal passage.
 - Lubricate the airway with the included jelly for easier insertion.
- **For Tension Pneumothorax:**
 - Only trained personnel should perform needle decompression using the ARS. Follow medical protocols for this procedure.

3. After Use:

- Dispose of used items according to local regulations.
- Replenish any used supplies to maintain the kit's readiness.

Disposal Instructions

- Dispose of all used medical supplies in accordance with local regulations.
- Follow guidelines for hazardous waste disposal for any items that may be contaminated.
- Ensure that the kit is restocked with fresh supplies after use.

Contact Information for Further Support

For any inquiries regarding the MICRO TRAUMA KIT NOW! ADVANCED SUPPLIES, please contact the manufacturer or refer to the product packaging for additional support information.

By following these safety instructions, you can ensure the effective use of the MICRO TRAUMA KIT NOW! ADVANCED SUPPLIES while adhering to safety guidelines. Always prioritize safety and seek professional medical assistance when necessary.

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