MICRO TRAUMA KIT NOW! BASIC SUPPLIES - MICRO TRAUM KIT NOW! BASIC SUPPLIES

The basic Micro Trauma Kit refill is not just another small tactical first aid kit. Each component was specifically chosen for its proven effectiveness in the field combined with never before seen packaging techniques to make the kit unobtrusive and completely packable for essential everyday carry. This trauma kit contains everything necessary to prevent deaths due to traumatic injuries. The Basic Kit contains the following: 1x Hemostatic dressing for wound packing/clotting 1x 4" Emergency Trauma Dressing) 6x 9" Medical Grade Easy Tape 1x Tourni-Kwik Compression Tourniquet 1x Large Heavy Duty Medical Gloves

Attributes

- Name: MICRO TRAUM KIT NOW! BASIC SUPPLIES
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100044831
- Mfr. No.: P-TKR-06
- Delivery weight: 0.195kg
- UPC: 812114025899

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MICRO TRAUMA KIT NOW! BASIC SUPPLIES Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! BASIC SUPPLIES. This kit is designed to provide essential first aid supplies for emergency situations. To ensure the safe and effective use of this product, please read and follow the safety instructions outlined in this guide.

General Safety Guidelines

- Always ensure the kit is stored in a cool, dry place, away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items.
- Keep the kit out of reach of children unless supervised by an adult.
- Use the kit only for its intended purpose: to provide first aid in emergency situations.
- Familiarize yourself with the contents of the kit and their specific uses before an emergency occurs.

Specific Safety Precautions for Use

- **Hemostatic Dressing:**
 - Use only on wounds that require packing or clotting.
 - Do not apply to nonbleeding areas.
- **Emergency Trauma Dressing:**
 - Ensure the dressing is applied tightly but not so tight as to cut off circulation.
 - Monitor the patient for signs of shock or distress.
- **Medical Grade Easy Tape:**
 - Use to secure dressings in place.
 - Avoid using on sensitive skin or open wounds.
- **TourniKwik Compression Tourniquet:**
 - Apply only if there is severe bleeding that cannot be controlled by direct pressure.
 - Follow the instructions on the tourniquet for proper application.
- **Medical Gloves:**
 - Always wear gloves when treating wounds to prevent infection.
 - Dispose of gloves properly after use.

Instructions for Installation and Usage

1. Preparation:

- Ensure the area is safe and the patient is calm.
- Gather all necessary supplies from the kit.

2. Using the Hemostatic Dressing:

- Remove the dressing from its packaging.
- Pack the dressing into the wound firmly and securely.
- Apply direct pressure to help control bleeding.

3. Applying the Emergency Trauma Dressing:

- Place the dressing over the wound.
- Secure it in place using the attached ties.
- Ensure that the dressing is tight but not restrictive.

4. Using the Medical Grade Easy Tape:

- Cut the required length of tape.
- Apply the tape over the dressing to hold it securely.

5. Applying the TourniKwik Compression Tourniquet:

- Place the tourniquet above the injury site (closer to the heart).
- Pull the strap through the buckle and tighten until the bleeding stops.
- Note the time of application.

6. Wearing Medical Gloves:

- Open the glove packaging without touching the inside.
- Put on the gloves before treating any wounds.
- Remove gloves carefully to avoid contact with contaminants.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not reuse any single-use items.
- Ensure that all waste is placed in a sealed container to prevent contamination.

Contact Information for Further Support

For any inquiries or further support regarding the MICRO TRAUMA KIT NOW! BASIC SUPPLIES, please refer to the manufacturer's contact details provided on the packaging.

This safety guide is provided to ensure the safe and effective use of your MICRO TRAUMA KIT NOW! BASIC SUPPLIES. For any safety concerns or to report unsafe products, please check the EU's Safety Gate platform for the latest updates. Thank you for prioritizing safety and preparedness.

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