# MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MULTICAM

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)

## Attributes

- Name: MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MULTICAM
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050369
- Mfr. No.: HW-TKN-MTKN-ESS-MC
- Size: 1 Person
- Delivery weight: 0.322kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650022

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# Micro Trauma Kit NOW! Essentials Supplies Safety Instructions

## Introduction

Thank you for choosing the Micro Trauma Kit NOW! (Micro TKN). This kit is designed to provide essential lifesaving medical supplies for everyday carry. It is crucial to understand how to use this product safely and effectively. This guide will help you familiarize yourself with the safety guidelines and proper usage of the Micro TKN.

## General Safety Guidelines

- Ensure that the Micro TKN is stored in a cool, dry place away from direct sunlight.
- Always inspect the kit before use to ensure all components are present and in good condition.
- Keep the kit out of reach of children and vulnerable individuals unless directly supervised by an adult.
- Do not use any expired medical supplies. Regularly check expiration dates and replace as necessary.
- Report any unsafe products or accidents to the appropriate authorities.
- Stay informed about recall updates on the EU's Safety Gate platform.

## Specific Safety Precautions for Use

- Trauma Gloves:
  - Use gloves to protect yourself when providing first aid.
  - Dispose of gloves properly after use to avoid contamination.
- Quickclot Bleeding Control Dressing:
  - Apply directly to the wound and hold in place until bleeding stops.
  - Do not use on nonbleeding wounds.
- Mini Compression Dressing:
  - Use for controlling bleeding and providing pressure.
  - Ensure the dressing is secure but not too tight to cut off circulation.
- Hypothermia Blanket:
  - Use the blanket to maintain body temperature in emergency situations.
  - Ensure the blanket is wrapped securely around the individual.
- Flat Fold Tape:
  - Use tape to secure dressings in place.
  - Avoid applying tape directly to the skin to prevent irritation.

## Instructions for Installation and Usage

### **1. Accessing Supplies:**

- The Micro TKN features pull tabs (BLIP) for easy access to supplies.
- Use one hand or a single finger to pull the tabs from either side of the kit.

### **2. Using the Kit:**

- In case of an emergency, retrieve the necessary supplies quickly.
- Follow the specific instructions for each item as outlined in the "Specific Safety Precautions for Use" section.

### **3. Storing the Kit:**

- After use, restock the kit with necessary supplies.
- Store the kit in its designated location for easy access.

## **Disposal Instructions**

- Dispose of used medical supplies according to local regulations.
- Do not dispose of medical waste in regular trash.
- Check with local health authorities for proper disposal methods for hazardous materials.

## **Contact Information for Further Support**

For any inquiries or concerns regarding the Micro Trauma Kit NOW!, please refer to the EUbased contact point provided with your product packaging. This ensures you receive accurate and timely assistance related to safety and usage.

Thank you for prioritizing safety with the Micro Trauma Kit NOW! Your preparedness can make a significant difference in emergency situations.

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