MICRO TRAUMA KIT NOW! SUPPLY REFILL - MICRO TRAUMA KIT ESSENTIAL SUPPLIES REFILL

The Micro Trauma Kit refill is not just another small tactical first aid kit. Each component was specifically chosen for its combat proven effectiveness combined with never before seen packaging techniques to make them as small as possible. The most popular MTKN configuration in service with police, military and medical providers worldwide. The ADVANCED kit is designed to treat traumatic injuries as per the MARCH protocol for those providers trained in TCCC concepts. The lightest and smallest fulltreatment trauma kit available. The Micro Trauma Kit supplies were designed from the ground up around a specific need for an ultra-concealable trauma kit. The Pro Version was designed primarily for First Responders, equipped to treat life threatening bleeding and bridge the time gap between injury and arrival of advanced medical providers. Paired with a tourniquet, the PRO kit allows for treatment of extremity and junctional bleeding with a minimal amount of training. The ESSENTIALS kit has just what you need to address the #1 cause of traumatic death - BLEEDING.

Attributes

- Name: MICRO TRAUMA KIT ESSENTIAL SUPPLIES REFILL
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100758
- Mfr. No.: TKR-MTKN-ESS
- Size: 1 Person
- Delivery weight: 0.177kg
- Shipping height: 38mm
- Shipping width: 171mm
- Shipping length: 197mm
- UPC: 810073650794

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Micro Trauma Kit Now! Supply Refill Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit Now! Supply Refill. This kit is designed to provide essential supplies for treating traumatic injuries effectively. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the kit is kept in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure all items are intact and within their expiration dates.
- Keep the kit out of reach of children and pets to prevent accidental use or ingestion of components.
- Always follow the instructions provided for each item in the kit to ensure proper usage.
- In the event of an emergency, call for professional medical help immediately.

Specific Safety Precautions for Use

- **Tourniquet Usage:** Ensure you are trained in the proper application of a tourniquet. Incorrect use can lead to further injury.
- **Wound Care:** Clean wounds thoroughly before applying any dressings to prevent infection.
- **Bleeding Control:** If treating severe bleeding, apply direct pressure to the wound and elevate the injured area if possible.
- **Training:** It is recommended that users receive training in Tactical Combat Casualty Care (TCCC) to maximize the effectiveness of the kit.
- **Age Restrictions:** This kit is intended for use by adults and trained individuals. It is not suitable for children without adult supervision.

Instructions for Installation and Usage

1. **Opening the Kit:** Carefully open the kit using the designated opening mechanism. Avoid using sharp objects that may damage the contents.
2. **Identifying Components:** Familiarize yourself with each component in the kit. A list of items and their uses is typically included in the kit.
3. **Using the Supplies:**
 - For bleeding control, use the appropriate dressings and tourniquet as per your training.
 - For minor injuries, use antiseptic wipes followed by adhesive bandages.
 - Follow the MARCH protocol for treating traumatic injuries:
 - M Massive Hemorrhage
 - A Airway
 - R Respiratory
 - C Circulation
 - H Hypothermia
4. **PostUse Care:** After using any supplies, ensure to dispose of them properly as per local regulations.

Disposal Instructions

- Dispose of used items in accordance with local waste disposal regulations.
- Do not dispose of medical supplies in regular household waste if they are contaminated.
- Consider using a sharps container for any needles or sharp objects included in the kit.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit Now! Supply Refill, please consult the product packaging for contact details of the manufacturer. Additionally, you can check for any safety updates or recall information on the EU's Safety Gate platform.

By following these safety instructions, you can ensure the effective and safe use of your Micro Trauma Kit Now! Supply Refill. Always prioritize safety and seek professional medical assistance when necessary.

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