MICRO TRAUMA KIT NOW! SUPPLY REFILL - MICRO TRAUMA KIT PRO SUPPLIES REFILL

The Micro Trauma Kit refill is not just another small tactical first aid kit. Each component was specifically chosen for its combat proven effectiveness combined with never before seen packaging techniques to make them as small as possible. The most popular MTKN configuration in service with police, military and medical providers worldwide. The ADVANCED kit is designed to treat traumatic injuries as per the MARCH protocol for those providers trained in TCCC concepts. The lightest and smallest fulltreatment trauma kit available. The Micro Trauma Kit supplies were designed from the ground up around a specific need for an ultra-concealable trauma kit. The Pro Version was designed primarily for First Responders, equipped to treat life threatening bleeding and bridge the time gap between injury and arrival of advanced medical providers. Paired with a tourniquet, the PRO kit allows for treatment of extremity and junctional bleeding with a minimal amount of training. The ESSENTIALS kit has just what you need to address the #1 cause of traumatic death - BLEEDING.

Attributes

- Name: MICRO TRAUMA KIT PRO SUPPLIES REFILL
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100759
- Mfr. No.: TKR-MTKN-PRO
- Size: 1 Person
- Delivery weight: 0.249kg
- Shipping height: 32mm
- Shipping width: 171mm
- Shipping length: 146mm
- UPC: 810073650800

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Micro Trauma Kit Now! Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit Now! This kit is designed to provide essential first aid supplies for treating traumatic injuries. It is crucial to understand how to use this kit safely and effectively. This guide outlines important safety information, usage instructions, and disposal guidelines to ensure your safety and the safety of others.

General Safety Guidelines

- Always read and understand the instructions before using the Micro Trauma Kit.
- Keep the kit out of reach of children to prevent accidental use.
- Regularly check the kit to ensure all components are present and in good condition.
- Store the kit in a cool, dry place, away from direct sunlight and extreme temperatures.
- If any component is damaged or expired, replace it immediately.
- Report any unsafe products or accidents to the relevant authorities.
- Check for recall updates on the EU's Safety Gate platform.

Specific Safety Precautions for Use

- **Hazard Identification:**
 - The kit contains medical supplies that may cause allergic reactions in some individuals.
 - Use caution when handling sharp objects, such as scissors or needles.
- **Avoiding Hazards:**
 - Do not use the kit for purposes other than its intended use.
 - Ensure that you are trained in basic first aid techniques before attempting to use the kit.
 - In case of severe injuries, always seek professional medical help immediately.
- **Age Specific Warnings:**
 - This kit is not intended for use by children without adult supervision.
 - Ensure that only trained adults handle the kit in emergency situations.

Instructions for Installation and Usage

1. Familiarize Yourself with the Kit:

- Open the Micro Trauma Kit and take note of all components included.
- Ensure you understand the purpose of each item in the kit.

2. Using the Kit:

- In case of a traumatic injury, assess the situation and ensure it is safe to approach the injured person.
- Follow the MARCH protocol for treating traumatic injuries:
 - **Massive hemorrhage:** Apply a tourniquet if there is severe bleeding.
 - **Airway:** Ensure the airway is clear.
 - **Respiration:** Check for breathing and provide assistance if necessary.
 - **Circulation:** Monitor blood circulation and treat for shock.
 - **Hypothermia:** Keep the person warm.
- Use the provided supplies as needed, following the instructions for each item.

3. PostUse:

- After use, restock any used items as soon as possible to maintain the kit's readiness.
- Clean any reusable items according to the manufacturer's instructions.

Disposal Instructions

- Dispose of used or expired medical supplies according to local regulations.
- Do not throw away medical waste in regular trash. Follow your local guidelines for hazardous waste disposal.
- Ensure that any sharp objects, such as needles, are disposed of in a designated sharps container.

Contact Information for Further Support

For any safety inquiries or further support regarding the Micro Trauma Kit Now!, please refer to the product packaging or manufacturer's website for contact details.

By following these guidelines, you can ensure safe and effective use of the Micro Trauma Kit Now! Always prioritize safety and seek professional medical help when necessary. Thank you for choosing our product.

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