MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-BLACK

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-BLACK
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100986
- Mfr. No.: TKN-NANO-ESS-BK
- Size: 1 Person
- Delivery weight: 0.147kg
- Shipping height: 44mm
- Shipping width: 83mm
- Shipping length: 133mm
- UPC: 810073652811

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Micro Trauma Kit NOW! Safety Instruction Guide

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Nano. This compact first aid kit is designed to provide essential medical supplies for emergency situations, particularly for controlling bleeding. Please read this safety instruction guide carefully to ensure safe and effective use of the kit.

General Safety Guidelines

- Ensure that the Micro Trauma Kit NOW! Nano is stored in a cool, dry place away from direct sunlight.
- Regularly check the kit for expired items or damage. Replace any items that are expired or compromised.
- Keep the kit out of reach of children unless under adult supervision.
- In case of an emergency, always call for professional medical assistance after administering first aid.
- Familiarize yourself with the contents of the kit and their uses before an emergency occurs.

Specific Safety Precautions for Use

- **Trauma Gloves:**
 - Use gloves to protect yourself and others from infection when treating wounds.
 - Dispose of used gloves properly and wash your hands immediately after use.
- **QuikClot Bleeding Control Dressing:**
 - Apply directly to the bleeding site and apply firm pressure.
 - Do not remove the dressing once applied, as it may disrupt clotting.
- **Pressure Device:**
 - Ensure the device is applied correctly to maintain consistent pressure on the wound.
 - Monitor the patient for signs of shock or distress while using the pressure device.

Instructions for Installation and Usage

1. Opening the Kit:

- Open the pouch carefully to avoid damaging any items inside.

2. Using Trauma Gloves:

- Roll out one pair of gloves from the kit.
- Put on the gloves before handling any wounds or applying medical supplies.

3. Applying QuikClot Dressing:

- Identify the bleeding source and clean the area if possible.
- Open the QuikClot dressing package and place it directly on the wound.
- Apply firm pressure for at least 3 minutes to control bleeding.

4. Using the Pressure Device:

- Position the pressure device over the dressing.
- Secure the device according to the instructions provided with it, ensuring it is snug but not cutting off circulation.

5. PostApplication:

- Once bleeding is controlled, seek professional medical help.
- Remove gloves and dispose of them properly. Wash hands thoroughly.

Disposal Instructions

- Dispose of all used medical supplies, including gloves and dressings, in accordance with local regulations.
- Do not dispose of medical items in regular household waste; check for local hazardous waste disposal options.
- Replenish the kit with new supplies as needed to maintain readiness for future emergencies.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW! Nano, please contact the manufacturer or your local distributor. Ensure you have your product details handy for a more efficient response.

Conclusion

The Micro Trauma Kit NOW! Nano is a critical tool for managing trauma situations effectively. By following these safety guidelines and instructions, you can ensure that you are prepared to respond to emergencies safely and efficiently. Remember, your safety and the safety of others is the top priority.

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