

# TRAUMA PAK W/QUICK CLOT - ADVENTURE MEDICAL KITS TRAUMA PAK W/QUIK CLOT

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

When it comes to emergencies, professionals know excess equipment only gets in the way. Act fast to provide immediate trauma care with this Trauma Pak which contains key supplies to control bleeding at the scene until medical professionals arrive. This essentials-only trauma pack contains Adventure® Medical Kit's QuikClot, a hemostatic dressing that stops bleeding three times faster by accelerating your body's natural clotting process with zeolite. In addition to the hemostatic dressing, this pack includes a combine pad, triangular bandage, and gauze dressings to help control the bleeding and duct tape to hold it all in place. Packed in a tough, waterproof pouch that easily fits into a BDU pocket or glove box, the Trauma Pak is compact, easy-to-use addition to your medical kit, first responder kit, or everyday carry bag that could help you save a life. Designed for Immediate Response Respond quickly to control bleeding until professional help arrives. This trauma kit contains the essentials you need to treat trauma in a quick-to-deploy format. Inside, you'll find a QuikClot® Advanced Clotting Sponge, trauma pad, triangular bandage, sterile gauzes, and other key supplies. Zeolite to Stop Bleeding Fast Control bleeding the hemostatic sponge that works on contact to accelerate your body's natural clotting process and stop bleeding three times faster. The sponge contains the hemostatic agent zeolite, an inorganic mineral that works quickly to induce clotting, in a non-adherent mesh bag. Critical Trauma Instructions Prepare yourself with the first aid information you need. This one-page guide includes life-saving trauma care instructions, including how to apply a tourniquet, stabilize fractures, manage gunshot wounds, and treat sucking chest wounds. Compact & Lightweight Have the tools you need to stop severe bleeding wherever you go. The Trauma Pak weighs only 4.4 oz. and is small enough to fit into a BDU pocket. Keep one in your glove box, range bag, backpack, first responder kit, and everyday carry bag so you're prepared when trauma occurs. Waterproof DryFlex™ Bag Have confidence this Trauma Pak will keep its vital contents ready to use at a moment's notice. The re-sealable DryFlex™ bag is completely waterproof and incredibly durable and can also be used to dispose of bio-waste or treat sucking chest wounds.

SUPPLY LIST

Bandage Materials	1 - Bandage, Conforming Gauze, 3"
	1 - Dressing, Gauze, Sterile, 2" x 2", Pkg./2
	1 - Dressing, Gauze, Sterile, 4" x 4", Pkg./2
Bleeding	1 - Gloves, Nitrile (Pair), Hand Wipe
	1 - 25 g QuikClot
	1 - Trauma Pad, 5" x 9"
Duct Tape	1 - Duct Tape, 2" x 26"
Fracture / Sprain	1 - Bandage, Triangular Wound Care
	4 - Antiseptic Wipe
Medical Information	1 - Trauma and Accident Management Instructions



## Attributes

- Name: [ADVENTURE MEDICAL KITS TRAUMA PAK W/QUIK CLOT](#)
- Manufacturer: [ADVENTURE MEDICAL KITS](#)
- Product no.: 100009128
- Mfr. No.: 2064-0292
- Size: 1 Person
- Delivery weight: 0.136kg
- Shipping height: 53mm
- Shipping width: 142mm

- Shipping length: 183mm
- UPC: 707708202926

## **Item details**

Made in USA

# Table of Contents

- [Startpage](#)
- [TRAUMA PAK W/QUICK CLOT Safety Instruction Guide](#)
- [About Us](#)

# TRAUMA PAK W/QUICK CLOT Safety Instruction Guide

## Introduction

Thank you for choosing the TRAUMA PAK W/QUICK CLOT by Adventure Medical Kits. This trauma pack is designed to provide essential supplies for controlling bleeding during emergencies. It is important to understand how to use this product safely and effectively. This guide will provide you with essential safety instructions and information to ensure safe use.

## General Safety Guidelines

- Ensure that the TRAUMA PAK is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the pack to ensure that all items are intact and within their expiration dates.
- Keep the trauma pack out of reach of children to prevent accidental use or ingestion of contents.
- Use the trauma pack only for its intended purpose, which is to provide immediate trauma care in emergency situations.
- Be aware of the potential hazards associated with using the items in the trauma pack, including allergic reactions to materials or improper application of dressings.

## Specific Safety Precautions for Use

- Before using any of the medical supplies, ensure that your hands are clean. Use the included hand wipes if necessary.
- If you are not trained in first aid, seek assistance from someone who is or call emergency services immediately.
- When applying the QuikClot hemostatic dressing, do not apply it directly to a wound without first ensuring that the area is clear of debris.
- In case of severe bleeding, apply direct pressure to the wound with the trauma pad before using other supplies.
- Follow the critical trauma instructions included in the pack carefully to ensure proper treatment of injuries.

## Instructions for Installation and Usage

### 1. Preparation:

- Open the waterproof DryFlex bag and remove the contents.
- Familiarize yourself with the items included in the pack: QuikClot sponge, trauma pad, gauze dressings, triangular bandage, duct tape, and antiseptic wipes.

### 2. Control Bleeding:

- For severe bleeding, apply the QuikClot sponge directly to the wound and hold it in place for at least 3 minutes to allow for clotting.
- If bleeding continues, apply a trauma pad over the QuikClot and secure it with duct tape.

### 3. Wound Dressing:

- Use the sterile gauze dressings to cover any additional wounds or to provide further protection.
- If necessary, use the triangular bandage to create a sling or secure dressings in place.

#### **4. FollowUp Care:**

- Once the bleeding is controlled, seek professional medical assistance as soon as possible.
- Monitor the injured person for signs of shock or other complications.

## **Disposal Instructions**

- Dispose of any used medical supplies in accordance with local regulations for biohazard waste.
- The DryFlex bag can be reused for storing additional supplies or for safe disposal of biowaste.

## **Contact Information for Further Support**

For any inquiries regarding the TRAUMA PAK W/QUICK CLOT, please refer to the packaging for the manufacturer's contact information or consult your local retailer.

## **Conclusion**

The TRAUMA PAK W/QUICK CLOT is a vital tool for emergency situations. By following these safety instructions and guidelines, you can help ensure that you and others are prepared to respond effectively to trauma. Always prioritize safety and seek professional medical help when necessary.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)