

# RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



## Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/TEK-LOK, COYOTE](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043722
- Mfr. No.: E10-3011B-CYT
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668267

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [RIGID TQ Case for SOF®TT/SOF®TTW Safety Instruction Guide](#)
- [About Us](#)

# **RIGID TQ Case for SOF®TT/SOF®TTW Safety Instruction Guide**

## **Introduction**

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This product is designed to provide secure storage and rapid access to your tourniquet in emergency situations. This safety instruction guide outlines essential safety guidelines, precautions, and instructions for proper use and maintenance of the RIGID TQ Case.

## **General Safety Guidelines**

- Ensure the RIGID TQ Case is used solely for its intended purpose, which is to hold the SOF®TT or SOF®TTW tourniquets.
- Regularly inspect the case for any signs of damage or wear. Do not use if damaged.
- Keep the case away from extreme temperatures and hazardous substances that may compromise its integrity.
- Ensure that the case is mounted securely to prevent accidental detachment during use.
- Always follow the manufacturer's guidelines for mounting and usage.

## **Specific Safety Precautions for Use**

- Ensure the tourniquet is properly secured within the case before use.
- When accessing the tourniquet, do so with one hand to maintain control of the case.
- Avoid using excessive force when pulling the tourniquet from the case to prevent injury.
- Be aware of your surroundings when deploying the tourniquet to avoid accidents.
- Ensure that the case is mounted at a height that allows for easy access without obstructing movement.

## **Instructions for Installation and Usage**

### **1. Mounting the Case**

- Choose a suitable location for mounting, such as a duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- For Police Duty Belt:
  - Use the BladeTech TEKLOK Belt Attachment, which fits up to a 2.25" belt.
- For MOLLE/PALS Attachment:
  - Use the Tactical Tailor Malice Clip for secure attachment.
- Ensure the case is mounted securely and at a height that allows for easy access.

### **2. Using the Case**

- To access the tourniquet, locate the exposed portion of the tourniquet.
- With one swift pull, release the tourniquet from the case.
- Practice using the case in a safe environment to ensure familiarity with its operation.

### **3. Securing the Hardware**

- Utilize the included Vibratite to secure the mounting hardware in place after achieving the desired configuration.
- Follow the instructions included with the Vibratite for proper application.

## **Disposal Instructions**

- Dispose of the RIGID TQ Case in accordance with local regulations regarding plastic and electronic waste.
- If the case is damaged or no longer needed, ensure it is disposed of in a manner that prevents environmental harm.

## **Contact Information for Further Support**

For any inquiries or concerns regarding the RIGID TQ Case, please refer to the manufacturer's contact information provided with your purchase. Ensure you have the product details ready for a more efficient response.

By following these safety guidelines and instructions, you can ensure the effective and safe use of the RIGID TQ Case for your SOF®TT/SOF®TTW tourniquets. Always prioritize safety and preparedness in emergency situations.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)