

RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, BLACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Supplies](#)

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.



Attributes

- Name: [ELEVEN 10 RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, BLACK](#)
- Manufacturer: [ELEVEN 10](#)
- Product no.: 100043725
- Mfr. No.: E10-3011M-BLK
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668243

Item details

Made in USA

Table of Contents

- Startpage
- RIGID TQ Case Safety Instruction Guide
- About Us

RIGID TQ Case Safety Instruction Guide

Introduction

Thank you for choosing the RIGID TQ Case for SOF®TT/SOF®TTW. This guide provides essential safety instructions and usage information to ensure you can effectively and safely utilize the product. Please read this guide thoroughly before use.

General Safety Guidelines

- Ensure that the RIGID TQ Case is used only for its intended purpose, which is to securely hold and allow rapid access to the SOF®TT or SOF®TTWide tourniquets.
- Regularly inspect the case for any signs of wear or damage. Do not use the case if it is compromised.
- Keep the case out of reach of children and ensure it is stored safely when not in use.
- Follow all local laws and regulations regarding the use of tourniquets and first aid supplies.

Specific Safety Precautions for Use

- **Friction Hold:** The case utilizes friction to hold the tourniquet securely. Ensure that the tourniquet is fully inserted and the case is closed before use.
- **OneHanded Deployment:** Practice the onehanded deployment of the tourniquet in a safe environment to familiarize yourself with the mechanism. This is crucial for effective use in emergency situations.
- **Blood and Wet Conditions:** Be aware that the case is designed for use even when hands are bloody or wet. However, ensure that you have a clear understanding of how to operate the case under these conditions.
- **Mounting Position:** When mounting the case, ensure it is securely attached to your duty belt, vest, or holster. Doublecheck the mounting hardware to avoid accidental detachment during use.
- **Avoid Modifications:** Do not alter or modify the RIGID TQ Case in any way that might compromise its functionality or safety.

Instructions for Installation and Usage

1. Mounting the Case:

- Choose a suitable mounting location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
- For Police Duty Belt: Use the BladeTech TEKLOK Belt Attachment, which fits belts up to 2.25 inches.
- For MOLLE/PALS Attachment: Use the Tactical Tailor Malice Clip.
- Ensure that the mounting hardware is securely fastened. Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained.

2. Inserting the Tourniquet:

- Open the RIGID TQ Case and insert the SOF®TT or SOF®TTW tourniquet into the case.
- Ensure that the tourniquet is positioned correctly and that the case is closed securely.

3. Deploying the Tourniquet:

- In an emergency, grasp the exposed portion of the tourniquet firmly.
- Pull swiftly to release the tourniquet from the case.
- Follow standard first aid protocols for applying the tourniquet once it is deployed.

Disposal Instructions

- Dispose of the RIGID TQ Case responsibly according to local regulations.
- If the case is damaged or no longer needed, ensure it is disposed of in a manner that does not pose a risk to others.
- Recycle materials where possible.

Contact Information for Further Support

For any questions or concerns regarding the RIGID TQ Case, please refer to the manufacturer's website or customer service resources. All products are designed, manufactured, and assembled in the USA, ensuring compliance with safety standards.

Conclusion

By following these safety instructions and guidelines, you can ensure the effective and safe use of the RIGID TQ Case for SOF®TT/SOF®TTW. Always prioritize safety and preparedness in emergency situations. Thank you for your commitment to responsible first aid practices.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk